



Product Information Sheet

Reference: _____

牌子/ Brand	YA-MAN	
型號 / Model	RF-BLOOM-HK	
原產地 /Place of manufacture	日本/JAPAN	
系統要求 / System Requirements	N/A	
主要功能 / Key features	<ul style="list-style-type: none"> ● RF 射頻輸出提升了 20%#，熱量更有效傳遞到皮膚角質層，激活膠原蛋白再生，達至提拉緊緻，淡化細紋效果。 ● 內置 Boost RF 程式，啟動後自動調節升溫，並感測皮膚溫度作最合適輸出，防止過熱灼傷。 ● 專利獨有設計雙環圓極射頻技術，配合極緻修護啫喱，讓射頻更加廣泛和均勻地傳遞到肌膚底層。 ● 機身輕巧細小，便於攜帶。 ● USB 充電設計，不受時間地域限制，隨時隨地都可以充電。 ● 建議護理時間: 每天只需 6 分鐘 (左、右臉各 3 分鐘)。 ● 每 3 分鐘會有提示聲 ● 啟動超過 6 分鐘會自動關閉電源。 ● 有效提升肌膚彈性 20%、輪廓提升 6 度、減淡法令紋 21%* (# 與 HRF-10T 比較) (*根據日本臨床實驗協會(JACTA)調查結果。實際使用效果因人而異，建議持續使用以達至更佳效果) <ul style="list-style-type: none"> ● RF output is increased by 20%#. With higher power, more heat can be reached to skin, resulting improved wrinkles and elasticity. ● The Boost RF program controls the efficiency of heat generation, prevent injury from overheating. ● Patent unique circular design warms broad skin area quickly ● Lightweight and compact design, easy for carrying. ● USB powererd design, easy to charge at anywhere, anytime. ● Recommended treatment period: 6 mins per day (Left and Right face 3 mins respectively) ● Sound notice by every 3 mins operation. ● Device will turn off automatically after 6 mins operation. ● Improve skin elasticity effectively by 20%, uplift contour by 6 degree, reduce smile lines by 21%* (# compared with HRF-10T) (*According to Japan Clinical Trail Association (JACTA) investigation results. The actual use effect varies from person to person, it is recommended to continue to use for better results.) 	
接頭 /Connectors	USB Type-C™	
螢幕 / Display	N/A	
配套 / Packaging List	<ul style="list-style-type: none"> ● AC 充電器 x 1 ● USB 線(USB Type-C™) x 1 ● 極緻修護啫喱 80g x 1 ● 使用手冊 x 1 <ul style="list-style-type: none"> ● AC Adapter x 1 ● USB Cable (USB Type-C™) x 1 ● Flawless Gel 80g x1 ● User Manual x 1 	
產品價格/ Price	請參照銷售備忘錄/發票 As per Sales Memo/Invoice	
附加費 / Additional Delivery Charge	N/A	
保修人身份 / Identity of repair service provider	Jebsen Consumer Products Service Centre 捷成消費品服務中心	
保修地址 / Repair service address	九龍旺角太子道西 193 號新世紀廣場第 2 座 19 樓 1917-1920 室 Unit 1917 – 1920, 19/F, Tower 2, Grand Century Place, 193 Prince Edward Road West, Mongkok, Kowloon	
換貨政策 / Exchange policy	你可於銷售單日期起計 14 天之內攜同流動通訊設備到我們的門市更換一次，詳情請參閱銷售單之條款及細則。 You may exchange the Mobile Equipment once within 14 days of the date of Sales Memo by returning the Mobile Equipment to our shops. Please refer to the terms and conditions on the Sales Memo for details.	
保修期 Warranty period	主機/ Main Unit	1 年 請登入公司網頁填妥保用登記表格: http://hk.jebsenconsumer.com/warranty-info-page For warranty registration, please visit and register online: http://hk.jebsenconsumer.com/warranty-info-page
	配件 / Accessory	N/A
詳細規格，請瀏覽 For details specification, please browse	http://hk.jebsenconsumer.com/hc-beaute-pure-hydrogen-water-bottle.html	

上述產品資料摘錄自生產商/供應商文獻及/或由生產商/供應商所提供之資料



The above information is extracted from literature and / or information on this product by manufacturer / supplier.
Hong Kong Telecommunications (HKT) Limited / CSL Mobile Limited 香港移動通訊有限公司
查詢熱線 Enquiry hotline: 2512 3123
<http://www.hkcsl.com>